


SINGLE USE TECHNOLOGY

by Rommelag



Healthcare Customer Event | Lyondellbasell Cologne, June 14+15, 2022 | Jakob Hansen & Stephan Gschwind



Table of contents

Overview: The Rommelag Group

- Part 1: Rommelag Flex Flecotec Flecozip, the Basic Idea Main Applications in the Pharmaceutical Industry Rommelag Single Use and Sustainability
- Part 2: Rommelag CMO Process Overview Rommelag CMO – Disposable filling system (DFS) Practical Use and Innovation

Summary / Key Message / Literature



ROMMELAG ENGINEERING


ROMMELAG CMO

Aseptic bottelpack filling systems and BFS inspection systems BFS contract manufacturing of liquid and semi-solid products **ROMMELAG** FLEX

Innovative containment solutions

ROMMELAG SERVICE

Individualised aftersales and Pharma Service






MAIN APPLICATIONS IN THE PHARMACEUTICAL INDUSTRY

Flexible Isolators for i.e. contained dispensing



Material Transfer System for contained bag charging & discharging



Vaccum Bags, for contained processes under pressure





ROMMELAG SINGLE USE AND SUSTAINABILITY


General Trend





SUT

General Pharma Trend



Do we have a conflict ?






ROMMELAG SINGLE USE AND SUSTAINABILITY

The Answer : LCA Study


LCA (Life Cycle Assessment): Cradle to grave study

ISO LCA Standard ISO 14040:2006 ISO 14044:2006





ROMMELAG SINGLE USE AND SUSTAINABILITY



LCA Study

Stainless Steel IBC




Single-Use FIBC





ROMMELAG SINGLE USE AND SUSTAINABILITY


Direct Impact Comparison:



Important Note: Datas valid for equipment used in the pharmaceutical (HP)API production /1/



Process Overview (cycle process) – What is BFS?



extrusion of parison with sterile air

and cutting of parison with sterile air overpressure

forming with sterile air or vacuum followed by filling

Container sealing



bottelpack[®] - how it works in general /2/3/4/5/



Link to BFS Process for Ampoules - YouTube



BOTTELPACK® WITH DISPOSABLE FILLING SYSTEM (DFS)


Conceptual design

- Optional dosing system to the **bottelpack**[®]
- Dosing of the DFS works with an peristaltic pump in place of time-pressure dosing
- Use of preassembled sterile disposable components in the product leading way, high sterility safeness
- Short reaction time, less time for construction activities by change of filling volume
- Design for small product batches from 100 ml to 5-10 litres and / or high-priced products (special in pre-clinic and clincial phases)
- Lower costs, lower barriers to entry for potential customers interested in a technical filling – quicker entry into the **bottelpack**[®] technology, more flexibility



BOTTELPACK® WITH DISPOSABLE FILLING SYSTEM (DFS)

Layout / overview main components



bag

Disposable adapter with sterile connector and connector for SIP

Manifold with preassembled disposable filling line

Pump unit, consisting of max. 12 peristaltic pumps on 312 – dimension 300 x 300 x 450 mm

Mandrel unit



BOTTELPACK® WITH DISPOSABLE FILLING SYSTEM (DFS)

Benefits

- Cleaning and cleaning validation costs (cross contamination) are greatly reduced when compared to stainless steel installations
- Lower cost, reduction in the amount of product required due to the elimination of sterile pipe work (high-priced products, small batches)
- Less format depending parts needed to cover a big width of different fill volumes
- Easy steps for conversion, high sterility safeness
- Independent of certain manufacturers of disposable parts works with different systems
- Ecological footprint related to a classic machine cleaning incl. cleaning validation is smaller
- Advantages in development = speed, few cavities, no cleaning effort good for freezing / thawing processes (biological drug products)



BOTTELPACK® WITH DISPOSABLE FILLING SYSTEM (DFS) PRACTICAL USE /6/7/8/9/

rhDNase (Pulmozyme, Genentech)

4 mg/ml formulation, 37°C for 15 min. Visual inspection, ELISA, CD, UV SEC, activity assay

Fully active, no aggregates, no permanent changes to conformational states 2-year refrigerated stability verified

EPO (Erytropoetin, Epoetin-Beta)

200 mycl with 500 bis 10`000 I.E. **BFS-Process with Purell 3020D** Visual inspection, ELISA, CD, UV SEC, activity assay

Fully active, no aggregates, no permanent changes to conformational states 2-year refrigerated stability verified

Attenuated Live-Virus Vaccines (flu vaccine and Rota-Virus vaccine)

0.2 and 2.3 ml **BFS-Process with Purell 1840H**

Recombinant RSV F Vaccine (Novovax)


No statistically significant differences in stability compared to conventional filled market products (glass & LDPE container) 2-year (flu vaccine) & 1-year (Rota-Virus) stability verified

The stability profile of RSV F vaccine under accelerated and intended storage conditions is comparable to the profile in glass vials and syringes



PRACTICAL USE AND POTENTIAL INNOVATION

Apiject



We bring together Blow-Fill-Seal and injection-molding technology so that pharmaceuticals and vaccines can be fill-finished in single-dose prefilled injectors at any scale.

ApiJect | Single-Dose, Prefilled Injectors using BFS Technology

UNITHER invests 68 million euros in the innovative project: Euroject® : Unither Pharmaceuticals (unither-pharma.com)

Both companies use blow-fill-seal technology as a foundation for their product development



SUMMARY / KEY MESSAGE

- Rommelag Flex's solutions ensure the highest personal and product protection when handling high potent APIs in the pharmaceutical industry.
- A disposable fill systems working with blow-fill-seal (BFS) technology would permit products fill (with limited volumes or small batch sizes) that could not normally use by BFS.



QUESTIONS?

I will use Google before asking dumb questions

I will use Google before asking dumb questions. www.mrburns.nl before asking dumb questions. I will use Google before asking dumb questions I will use Google before asking dumb questions. I will use Google asking dumb questions. I will use Google before asking dumb questions. I will use Google before asking dumb questions. I will use Google asking dumb questions. I will use Google before asking dumb questions. I will use Google before asking dumb questions. I will use Google asking dumb questions. I will use Google before asking dumb questions.

-

ROMMELAG CMO

LITERATURE / REFERENCE

- 1. Mohammed Gaffar Khan, Master thesis from Institute for Polymer Science and Processing (iPSP) Hochschule Aalen für Technik und Wirtschaft, comparing environmental aspects of traditional stainless-steel technology and single use plastic technology in API & pharmaceutical production and conception of API packaging, Jan. 2018
- 2. European Commission, EU Guidelines to Good Manufacturing Practice. Annex 1, Manufacture of Sterile Medicinal Products (Brussels, Nov. 2008).
- 3. FDA, Guidance for Industry. Sterile Drug Products Produced by Aseptic Processing Current Good Manufacturing Practice (Rockville, MD, Sept. 2004).
- 4. R. Oschmann, and O.E. Schubert, Eds, Blow-Fill-Seal Technology, (CRC Press, Stuttgart, 1999).
- 5. PDA, Technical Report No. TR77, The Manufacture of Sterile Pharmaceutical Products Using Blow-Fill-Seal Technology, (2017)
- 6. Verjans, B. Reed, C. (2012). "Assessing Filling Technologies for Contamination Risk." Biopharm International. 25(3), pp. 46-58.
- 7. Gregory Bleck, Ph.D. Global Head R&D Biologics Madison, Wisconsin, Biologics Case Study: Compatibility assessment of a model monoclonal antibody formulation with ADVASEPT ® and glass container system, 2015 Catalent Pharma Solutions
- 8. Darin Zehrung, Next-Generation Vaccine Delivery Technology Meeting Geneva, Switzerland, Feb. 2014
- 9. Poster Novovax Inc., Feasibility Evaluation of Blow Fill Seal Process and Compatibility with Aluminium Phosphate Adjuvanted Recombinant RSV F Vaccine



A



Jakob Hansen

Managing Director

Rommelag Flex
In der Eschenau 5
74405 Gaildorf
Jakob.Hansen@rommelag.com

in

+49 7971 256 0

Stephan Gschwind
Managing Director
Maropack AG
Industriestrasse Briseck 4
6144 Zell, Switzerland
Stephan.Gschwind@rommelag.com
₩ +41 41 989 74 00

- www.rommelag.com
- www.linkedin.com/company/rommelag
- www.youtube.com/c/RommelagGroupofCompanies